Excel Avensis

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User manual

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Inhoudsopgave	Inhaltsverzeichnis
Table of contents	Table des matières
Índice	

Nederlands.....	4
English.....	26
Deutsch.....	48
Français.....	70
Español.....	100

Nederlands

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English

Table of contents

General introduction	30
1. Identification	32
2. General information	33
2.1 Components and structure of your rollator	33
2.2 General explanation about your rollator	33
2.3 Intended use	33
2.3.1 Indication and contraindication	34
3. Safety instructions	34
3.1 General safety instructions	34
3.2 Warnings for safe use	34
4. Specifications	36
5. Components of your Excel Avensis rollator	37
5.1 Folding and unfolding the rollator	37
5.1.1 Folding the rollator	37
5.1.2 Unfolding the rollator	37
5.2 The brakes	37
5.2.1 Braking	38
5.2.2 Parking brake	38
5.2.3 Adjusting the braking force	38
5.3 The seat	39
5.3.1 Adjusting the backrest	39
5.4 The steppers	39
5.5 The push handles	40
5.6 The storage bag	40
5.7 The wheels	41

6.	User instructions	42
6.1	Balance.....	42
6.1.1	Walking with your rollator.....	42
6.1.2	Sitting	42
6.2	Navigating obstacles	43
6.2.1	Climbing an obstacle	43
6.2.2	Descending an obstacle	43
6.3	Navigating slopes.....	44
6.3.1	Ascending a slope.....	44
6.3.2	Descending a slope.....	44
6.3.3	Moving surfaces	44
7	Maintenance	45
7.1	The tyres.....	45
7.1.1	PU tyres	45
7.2	The brakes	45
7.3	The frame	45
7.4	The rear wheels	45
7.5	The swivel wheels and the front fork	46
7.6	The upholstery	46
7.7	Cleaning your rollator	46
7.8	Periodic maintenance of your rollator.....	47
7.9	Reuse.....	47
8.	Troubleshooting and solutions.....	48
9.	Warranty	49
9.1	Warranty application	49
9.2	Warranty conditions	49

General introduction

Before using your new Excel Avensis, you and any potential caregiver should thoroughly read and understand this user manual. By reading and understanding this manual, you also agree to the conditions mentioned in this document.

You are now the owner of a quality product from Excel Mobility. We thank you for placing your trust in our Excel® products. The Excel Avensis is truly a high-quality product. Excel Mobility is committed to continuously improving the quality and reliability of all our products. Therefore, we reserve the right to make changes to this user manual without prior notice.

It is important to carefully read this manual before using your rollator. This manual contains essential information about safely using and maintaining your rollator. We recommend keeping this manual safe, as it may be useful later.

The safety instructions in this manual are general guidelines intended to provide broad directions. You may develop your own methods for performing common tasks. However, we strongly advise consulting a professional to develop safe and effective techniques tailored to your daily activities and physical abilities.

Your new rollator requires regular and systematic maintenance. While much of this can be performed by yourself, we recommend having your rollator inspected by a professional at least once a year.

Note!

In this user manual, you will find tips and warnings clearly marked with symbols and specific text formatting as shown below:

Tip

Tips included in this manual are designed to help you make the best use of your Excel Avensis rollator.

Warning

These warnings must always be followed to avoid damage to the Excel Avensis rollator or injury to yourself.

Fill in the details of your authorised dealer and product below:

Company name:

Address:

Telephone number:

Email address:

Website:

Date of purchase:


Stamp or sticker from your authorised dealer (if applicable):

If, in addition to the frame label on the product, a second separate frame label was provided with your product, you can attach it below for easier reference later:

1. Identification

Your rollator is equipped with a frame label, located on the cross frame beneath the seat. The frame label contains important information about your rollator, such as the product name, manufacturer, authorised European representative, unique device identifier (UDI), and technical specifications. Below is an example of a frame label for reference. The actual label on your rollator may differ. The details on the label correspond specifically to your rollator.

Van Os Medical B.V. is the importer of this Excel Mobility product. The product is manufactured by Excel Mobility Manufacturing Co. Ltd and is represented in Europe by SUNGO Europe B.V.

1. Model name + description:

The model name of your rollator starts with the brand name “Excel” followed by an additional model description. For any inquiries or more information about your rollator, it is important to have the additional model description at hand, in this case, “Excel Avensis”.

2. Production date:

This indicates when your rollator was manufactured.

3. Serial number:

Each rollator has a unique serial number. You will need this number if you have technical questions or wish to order parts for your rollator.

Warning

Any structural modifications to the Excel rollator will void the warranty entirely.

Warning

Your Excel rollator should only be used for its intended purpose.

2. General information

2.1 Components and structure of your rollator

1. Reflector on the handle
2. Brake lever
3. Adjustable push handles
4. Fixation clip to secure the push handles
5. Adjustment knob for brake force
6. Rear wheel
7. Handle
8. Backrest
9. Seat
10. Storage bag
11. Front wheel

The components shown here are part of the standard configuration.

Photo 1

Familiarise yourself with these components before studying this user manual further. Designs and specifications are subject to change without prior notice.

2.2 General explanation about your rollator

Your new Avensis rollator is made of ultra-lightweight yet strong aluminium, making it exceptionally manoeuvrable and durable. The rollator features integrated safety functions such as built-in brake cables and automatic locking when folded. The rollator is suitable for indoor and outdoor use. With ergonomically adjustable handles, soft shock-absorbing wheels, and a comfortable backrest, you will experience maximum comfort both when walking and taking a break.

2.3 Intended use

The rollator is intended for individuals who require support while walking or who have reduced mobility. The user must be able to operate the rollator independently and safely. Do not use the rollator if physical or mental limitations make safe operation impossible. If in doubt, always consult a doctor or specialist before using the rollator.

Tip

Read this manual thoroughly for optimal use of your rollator.

2.3.1 Indication and contraindication

Indication:

- Neurological conditions: Parkinson's disease, multiple sclerosis, consequences of a stroke.
- Musculoskeletal conditions: osteoarthritis, arthritis, osteoporosis.
- Cardiovascular and respiratory conditions: heart failure, COPD.
- Balance disorders and dizziness.
- Postoperative conditions: after surgery on bones or joints, such as knee or hip surgery.
- Reduced muscle strength: weakness in the legs or hips.
- Walking insecurity: reduced gait stability or fear of falling.

Contraindication:

The rollator must not be used when:

- balance is impaired to such an extent that the user cannot operate the rollator safely (risk of falling).
- insufficient perceptive ability prevents proper assessment of the surroundings and the rollator.
- the user is unable to sit down on the rollator, which significantly limits its usability.
- the user is physically unable to operate or manoeuvre the rollator.

3. Safety instructions

The Excel Avenis rollator is an extremely safe and stable product when the guidelines outlined in this user manual are followed. However, improper use may result in dangerous situations. Van Os Medical B.V. accepts no liability for personal injury and/or property damage if the recommendations and warnings in this manual are not followed.

3.1 General safety instructions

Regularly check your rollator to ensure it is in good condition. Malfunctioning parts may pose safety risks. Periodic inspection, correct adjustments, and timely replacement of worn or broken components will enhance the product's lifespan. An authorised Excel dealer, using only Excel replacement parts, will ensure the rollator remains in optimal condition through preventive maintenance.

3.2. Warnings for safe use

To prevent damage to your property, the product itself, or personal injury, adhere to the following warnings. Van Os Medical B.V. is not liable in any form. By reading these warnings, you agree to the terms stated in this manual. Any modifications not described in this manual will void the warranty entirely.

- Always check upon receipt that all rollator components are present and undamaged. If any components are missing or damaged, contact your dealer immediately.
- Before using the rollator, always verify that the brakes are functioning properly. Squeeze both brake levers and attempt to push the rollator forward. If the rear wheels slide rather than roll, the brakes are working correctly.
- Use the parking brakes when stopping or parking the rollator.
- Always use an elevator to transport the rollator to a higher or lower level.
- When using the rollator in an elevator, it must always be in the braking position.
- Engage the parking brake before sitting.
- Ensure that all wheels remain in contact with the ground during use to maintain proper balance.
- Adjust the rollator to the correct height to prevent improper posture. Refer to section 5.5 for details.
- When unfolding the rollator, confirm that the frame's locking mechanism has clicked securely into place, see figure 3c
- Be cautious of clothing or body parts that may become trapped when folding or unfolding the rollator.
- The storage bag can hold a maximum of 5 kg. Exceeding this weight may alter the rollator's centre of gravity, increasing the risk of tipping.
- When modifications are made to the rollator that are not described in this user manual, the warranty rights are entirely voided.

- Do not use the rollator on motorways; use it on pavements or accessible roads whenever possible.
- Never use your rollator if it is not functioning properly or if any of its components are malfunctioning or missing.
- Do not use your rollator if the brakes are not working.
- Do not stand on your rollator.
- Avoid using the rollator on uneven surfaces, such as sand, gravel, or wet and slippery terrains.
- Never ascend kerbs directly; always use ramps when available.
- Ensure the rollator is not overloaded beyond the stated maximum user weight of 136 kg, including bags or additional load.
- Do not use the rollator for transporting heavy objects or people.
- Do not move the rollator while seated or allow others to do so.
- Keep small children from playing with the rollator.
- Never attach objects to the wheels, as this could cause tipping or damage, compromising safety.
- Avoid using the rollator on escalators or moving walkways; see section 6.3.3 for more information.
- Do not push the rollator forward while it is in the braking position. This may damage the wheels and pose a safety hazard.
- Do not push the rollator too far ahead of you, as this may cause falls or improper positioning. See figure 3a.

Fig. 3a

Fig. 3b


Fig. 3c

4. Specifications

	S	M
Total length (cm)	67	67
Total width (cm)	58	58
Total height (cm)	80-90 (adjustable push handles)	80-95 (adjustable push handles)
Total weight (kg)	6	6,1
User weight (kg)	136	136
Seat width (cm)	45	45
Seat depth (cm)	20	20,5
Seat height (cm)	54	60
Maximum load of storage bag (kg)	5	5
Distance between handles (cm)	45,5	45,5
Turning radius (cm)	54,5	54,5
Recommended user height (cm)	150-165	165-175

	S	M
Length when folded (cm)	67	67
Width when folded (cm)	20	19
Height when folded (cm)	82	85

5. Components of your Excel Avenis rollator

In this chapter, the various components of your rollator are described. The rollator includes several components that are adjustable or can be folded. If you wish to attach or remove components, adjust or reposition them, always follow the instructions described in this user manual.

5.1 Folding and unfolding the rollator

If you want to fold or unfold the rollator, there are a number of steps to follow. Always exercise caution when folding or unfolding the rollator to avoid damage.

5.1.1 Folding the rollator

To fold the rollator, follow these steps:

- Stand behind the rollator and remove the storage bag (see section 5.6).
- Pull the strap attached to the seat upwards. This will automatically fold the rollator (see photos 2–5).

Photo 2

Photo 3

Photo 4

Photo 5

5.1.2 Unfolding the rollator

To unfold the rollator, follow these steps:

- Gently press the seat tubes downwards to unfold the rollator. You will hear a confirming click once the rollator is fully unfolded (see photos 6–8).
- Replace the storage bag, and the rollator is ready for use (see photo 9).

Photo 6

Photo 7

Photo 8

Photo 9

5.2 The brakes

The Excel Avenis is equipped with two types of brakes: a parking brake and a brake for deceleration. To use the brakes, follow these steps:

5.2.1 Braking

To slow down while walking behind the rollator, squeeze the brake lever (see photo 10). Avoid pushing the rollator forward while braking, as this may damage the tyres and brakes.

Photo 10

5.2.2 Parking brake

- Push the brake lever downwards to engage the parking brake (see photo 11). You will hear a click when it is locked.
- To release the parking brake, squeeze the brake lever. You will hear another click as the parking brake is disengaged.

Photo 11

5.2.3 Adjusting the braking force

An adjustment knob near the rear wheel allows you to adjust the braking force. Push the adjusting knob downwards (see photos 12-13) and turn it to the left to increase the braking force or to the right to decrease it (see photo 14).

Photo 12

Photo 13

Photo 14

5.3 The seat

The Excel Avensis is equipped with a seat for short periods of comfortable sitting. To use the seat, follow these steps:

- Ensure the rollator is in the parking brake position (see section 5.2.2).
- Once the parking brake is engaged, sit down on the seat.
- Ensure your feet are positioned between the rear wheels and that you are sitting with your back facing the front wheels (see photo 15).

Photo 15

5.3.1 Adjusting the backrest

The backrest can be adjusted in both depth and height.

To adjust the depth, press the silver button and pull the backrest to slide it to the desired position (see photos 16–17).

Photo 16

Photo 17

Photo 18

Photo 19

To adjust the height, turn the black star-shaped knob and move the backrest up or down as needed (see photos 18–19). Ensure the knob is securely tightened before sitting to prevent the backrest from slipping down.

5.4 The steppers

The Excel Avensis is equipped with two steppers, which can be used to apply additional leverage when navigating over a raised surface. However, caution is required when using the steppers.

To use the steppers, follow these steps:

- Place your foot on the stepper to lift the front wheels of the rollator (see photo 20).
- This creates a lever effect, causing the front wheels to tilt upwards.
- Navigate over the raised surface and remove your foot from the stepper to allow the front wheels to return to the ground.

Photo 20

5.5 The push handles

The push handles of the rollator are height-adjustable. To adjust the height of the push handles, follow these steps:

- Loosen the fixation clip (see photo 21).
- Press the adjustment buttons located at the back of the rollator frame (see photo 22).
- Pull the push handle upwards to the desired height (see photos 23–24). If the push handle is too high, press it downwards. Once the push handle is at the desired height, release the adjustment buttons, and they will automatically lock into place with an audible click.
- Re-secure the fixation clip (see photo 25).

Photo 21

Photo 22

Photo 23

Photo 24

Photo 25

5.6 The storage bag

The Excel Avensis is equipped with a storage bag. To remove the bag, follow these steps:

- First, detach the Velcro strap (see photo 26).
- Slide the loop off the rollator (see photos 27–28).
- To reattach the bag, slide the loop back onto the rollator and secure the Velcro strap.

Photo 26

Photo 27

Photo 28

Warning

The storage bag has a maximum weight capacity of 5 kg. Exceeding this limit may alter the rollator's center of gravity, increasing the risk of tipping.

5.7 The wheels

The rollator is fitted with front and rear wheels. The front wheels measure 9 inches, while the rear wheels are 7.5 inches. Regularly inspect the wheels for wear and tear.


Photo 29

Photo 30

The front wheels are attached to the frame via a front fork, which is critical for steering the rollator. If steering becomes difficult or the front wheels wobble, the wheels may need adjustment.

6. User instructions

It is important to learn the correct and safest way to use the rollator. Carefully study these user instructions and familiarise yourself with all functions and components of the product. Always consult a professional for advice on developing the proper and safest methods suited to your physical abilities.

6.1 Balance



Everyday activities, such as walking with or sitting on the rollator, affect its balance. This relates to the distribution of weight and the shifting of the center of gravity. Follow the instructions below to prevent the rollator from tipping or other dangerous situations.

This user manual describes guidelines for independent users without a caregiver. If you can no longer use the rollator independently, we recommend stopping its use to prevent accidents.

6.1.1 Walking with your rollator

When walking with your rollator, take the following points into account:

- Always maintain an upright posture. A hunched posture is detrimental to your back (see figures 4a and 4b).
- To achieve the correct posture, the push handles must be adjusted to the proper height. To find the correct height, let your arms hang naturally by your sides and adjust the handles to the height of your wrists (see figure 4c). Instructions on adjusting the push handles can be found in section 5.5

6.1.2 Sitting

When you want to sit on the rollator, for example, after a long walk, it is crucial to engage the parking brake. This prevents the rollator from moving and allows you to sit safely. Instructions on how to engage the parking brake can be found in section 5.2.2.

6.2 Navigating obstacles

You may encounter situations where you need to navigate over or around obstacles with your rollator. This manual provides instructions on how to do this safely.

6.2.1 Climbing an obstacle





When climbing an obstacle independently, do so in a forward motion.

- Press down on one of the steppers with your foot to lift the front wheels using the lever effect (see figures 4d and 4e).
- Carefully place the front wheels of the rollator on the obstacle, push the rollator forward, and follow it yourself (see figures 4f and 4g). This is the correct way to climb obstacles.
- Ensure that no items are placed on the seat when climbing an obstacle, as they may fall off (see figure 4h).

6.2.2 Descending an obstacle

When descending an obstacle independently, do so in a forward motion.

- Roll the front wheels onto the lower surface and apply light braking to carefully descend the obstacle with the rollator.
- Once both you and the rollator have fully descended the obstacle, release the brake and continue walking.
- Ensure that no items are placed on the seat when descending an obstacle, as they may fall off (see figure 4i).

 **Warning**
Never attempt to navigate obstacles higher than 5 cm.

6.3 Navigating slopes

It is possible to encounter situations where you need to ascend or descend slopes. When this occurs, we recommend always moving in a forward and smooth motion.

6.3.1 Ascending a slope

When ascending a slope independently, we advise leaning slightly forward. This shifts your center of gravity forward, making it easier to ascend the slope. Avoid placing additional weight, such as bags, on your rollator during this process, as this increases the risk of tipping (see figures 4h and 4i).

Always aim to ascend directly toward the highest point of the slope. Ascending diagonally increases the risk of falling.

6.3.2 Descending a slope

When descending a slope, always face forward. Never lean forward, as this may cause the rollator to tip. To increase stability, it is better to lean slightly backward. Always aim to descend in a straight line. If you need to stop while descending, use the rollator's parking brake. To maintain control, do not descend too quickly; regulate your speed by applying light braking when necessary.

Tip

Carefully assess a slope or obstacle before attempting to navigate it. Knowing what to expect allows you to proceed safely.

6.3.3 Moving surfaces

The Excel Avenis is not suitable for use on dynamic surfaces such as escalators or moving walkways. For safety reasons, we advise only using your Avenis on fixed, stable surfaces designed for this purpose. Always use an elevator to transition between floors or levels.

7 Maintenance

Your Excel Avensis rollator requires periodic maintenance. This is essential for a long lifespan and optimal user comfort. A poorly maintained rollator is more likely to develop technical problems, will function less smoothly, and will no longer be covered under warranty. Preventive maintenance is therefore very important. We recommend having your rollator checked by a qualified Excel Mobility dealer at least once a year. The dealer will use only Excel replacement parts if necessary. Regular inspections ensure that your rollator continues to perform perfectly for many years.

You can also do a lot yourself to keep your rollator in optimal condition. Regularly checking the rollator and performing minor maintenance when needed will extend its lifespan and improve your user experience. The following sections describe this minor maintenance in more detail.

7.1 The tyres

The tyres of your rollator should be checked regularly for wear. If the wear is advanced, the tyres need to be replaced.

7.1.1 PU tyres

Polyurethane (PU) tyres are puncture-proof. These tyres are identifiable by the absence of a valve. The advantage of PU tyres is that they do not need to be inflated and cannot puncture. However, PU tyres are still subject to wear. Therefore, check the tread weekly to ensure it is still sufficient.

7.2 The brakes

Inspect the brakes for sharp edges. If sharp edges are present, the brakes must be replaced. Also check that all brake components are securely fastened. Frequent use may cause parts to loosen. Additionally, ensure that the braking mechanism operates smoothly. If it does not, apply a few drops of WD40 oil or Teflon to the brake pivot points. Finally, make sure the brakes do not come into contact with the tyres while the rollator is in motion.

7.3 The frame

The frame of your rollator forms the foundation of the product and must therefore remain in good condition. Check the frame for the following points:

- Check whether the rollator folds and unfolds without resistance.
- Check whether the rollator moves forward in a straight line.
- Check whether the frame is bent or damaged.

All metal surfaces of the frame are coated with an electrophoretic paint layer that provides long-term protection against corrosion. The fastening materials, such as screws, are also corrosion-resistant. Regularly check whether the coating is intact and have any damage repaired to maintain corrosion protection.

If any of the above points show irregularities, contact the dealer where you purchased the rollator. Only your dealer is authorised to carry out frame repairs.

7.4 The rear wheels

You can also inspect the rear wheels yourself. Proper adjustment of the rear wheels enhances the rollator's stability and reduces the effort required to move it. Check the rear wheels for the following:

- Ensure the rear wheels do not rub.
- Check for play in the rear wheels. If present, it is often due to loose wheel housings.

If any of the above points show irregularities, contact the dealer where you purchased the rollator. Only your dealer is authorised to carry out repairs to the rear wheels.

7.5 The swivel wheels and the front fork

The swivel wheels at the front of the rollator and the fork to which they are attached must move smoothly for optimal steering. Inspect the swivel wheels and fork for the following:

- Ensure the nut is neither too loose nor too tight. A loose nut can cause the swivel wheels to wobble or come off, while a tight nut can make steering difficult. If the nut is too tight, the bearings will be compressed, preventing the balls from moving freely.
- Inspect the fork housing for signs of wear.

If any of the above points show irregularities, contact the dealer where you purchased the rollator. Only your dealer is authorised to carry out repairs to the swivel wheels or fork.

Tip

Check that the nut on the swivel wheel is properly fastened by spinning the wheel. It should come to a gradual stop.

7.6 The upholstery

The upholstery of your rollator is also an important component. Torn or worn upholstery can no longer support human weight, leading to dangerous situations. Regularly inspect the upholstery for holes, tears, or worn areas.

If the upholstery shows signs of damage, contact the dealer where you purchased the rollator. Only your dealer is authorised to carry out repairs to the upholstery.

7.7 Cleaning your rollator

Cleaning your rollator is important and should be done regularly. Follow these steps for cleaning:

- Use mild soap and water to clean the upholstery, chassis, and plastic components. Never use abrasive cleaners, as these can damage the finish.
- Never use a high-pressure washer or steam cleaner.
- Treat the chassis finish regularly with wax. However, do not use solvent- or abrasive-based wax, harmful chemicals, or silicone spray.
- Always dry the rollator thoroughly after cleaning or exposure to rain.

7.8 Periodic maintenance of your rollator

Preventive maintenance is very important, and many tasks can be done by yourself (or with assistance). Below is a checklist of maintenance intervals for the components described above.

Monthly

- Check whether the rollator folds and unfolds smoothly.
- Inspect the upholstery for damage or wear.
- Ensure there is no excessive play in the front and rear wheels.
- Test the brakes for proper functionality.

Quarterly

- Check that all bolts and nuts are securely fastened, and tighten if necessary.
- Inspect the tread on the tyres for wear.

7.9 Reuse

The rollator is suitable for reuse. Before being handed over to another user, the rollator must be inspected and hygienically cleaned by an authorised Excel Mobility dealer. When transferring the rollator, remember to provide this user manual and the dealer's contact details to the new user.

Before reuse, the following checks must be carried out by an authorised Excel Mobility dealer:

- Visual inspection for external damage
- Visual inspection for correct assembly of all components
- Functional check of all manual functions
- Complete inspection
- Cleaning and disinfection

Hygiene measures for reuse:

Before the rollator can be used again, it must be hygienically cleaned. All surfaces that come into contact with the user must be cleaned using a suitable disinfectant. Disinfectants must be diluted and applied according to the manufacturer's instructions. Always read the label before use.

8. Troubleshooting and solutions

During use, you may encounter problems with your Excel rollator. The possible issues are described in the table below, along with their solutions. If you experience problems with your Excel rollator that are not listed in the table, we recommend contacting an authorised Excel Mobility dealer.

Problem	Solutions			
	1	2	3	4
There is play in the rollator	X			
The rollator squeaks and/or rattles	X			
The swivel wheels wobble	X	X		X
The rollator is hard to move			X	
The rollator has slow steering	X	X		X
The rollator pulls to the left		X		X
The rollator pulls to the right		X		X

1	Check that all bolts and nuts are securely fastened.
2	Ensure the rear and swivel wheels are adjusted equally.
3	Check if the bearings need replacement.
4	The front forks of the swivel wheels may be misaligned. Check if they are properly adjusted.

9. Warranty

Your proof of purchase serves as your warranty certificate.

9.1 Warranty application

Your rollator comes with a manufacturer's warranty. However, your supplier may also offer an additional warranty. This manual exclusively covers the manufacturer's warranty provided by Excel Mobility.

The manufacturer's warranty is only granted to you as a consumer. It cannot be used for commercial purposes, such as rental or institutional use. The warranty is limited to defective materials and hidden defects. Excel Mobility guarantees the following warranty periods for your Excel Avensis rollator:

- Warranty period: side and cross frames, rims, bearings, brakes, forks, and all components not specifically mentioned below: 2 years.
- No warranty on: upholstered components and other wear-and-tear parts, such as the seat, backrest, tyres, paint, and coating.

9.2 Warranty conditions

The warranty period begins on the date of purchase of your rollator. To claim warranty on defects within the warranty period, you must present the rollator for repair to your authorised dealer, along with your proof of purchase.

Repairs and replacements are carried out by an authorised Excel Mobility dealer. If the rollator is used outside the area covered by the authorised dealer responsible for the warranty, repairs can be performed by any other dealer designated by the manufacturer.

If a part of the rollator covered by the warranty needs repair or replacement due to a manufacturing or material defect, and the rollator is owned by the original purchaser, the relevant parts will be repaired or replaced free of charge by the authorised dealer.

The warranty does not cover labour costs incurred during repair or replacement. Any replaced or repaired part is subject to the same warranty conditions as the rollator itself. Wear-and-tear parts are generally not covered under the warranty unless their wear is directly caused by a manufacturing defect.

Excel Mobility is not responsible for replacing or repairing the rollator if the following conditions are met:

- The rollator and its components are not maintained in accordance with the manufacturer's recommendations or original parts are not used.
- The rollator or its components are damaged due to negligence, accident, or improper use.
- The rollator or its components are modified from factory specifications, or repairs are conducted before notifying the service agent.
- The product does not have an original factory frame label (see Chapter 1).

The rollator described and shown in this user manual may differ slightly in detail from your own model. However, all instructions are relevant regardless of any differences in detail. We reserve the right to make changes to the weight, dimensions, or other technical data described in this user manual without prior notice. All drawings, dimensions, and capacities shown in this user manual are approximate and may not fully match the given specifications.

For warranty service, contact the authorised dealer where you purchased your rollator. If you are not satisfied with the assistance provided by your dealer, you can contact Excel Mobility. The importer's contact details can be found on the back of this manual.

Warning

Excel Mobility cannot be held liable for any consequential or incidental damage. Although this user manual has been carefully compiled, it is not comprehensive. If you intend to use your rollator in a manner not prescribed in this manual, we recommend discussing this with an authorised Excel Mobility dealer beforehand.

Warning

The warranty is only valid within the indicated period. Any modifications to the rollator that affect its structural integrity will result in complete loss of warranty. Contact Excel Mobility for detailed warranty and delivery terms, as well as a list of authorised dealers.

Tip

Maintenance is not covered under the warranty. Your dealer may deviate from the recommended maintenance intervals.

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Op dit moment maakt u kennis met een kwaliteitsproduct van Excel Mobility. Wij willen u bedanken voor het vertrouwen in onze Excel © producten. De Excel Avensis is een echt kwaliteitsproduct. Het beleid van Excel Mobility is erop gericht om doorlopend de kwaliteit en betrouwbaarheid van al onze producten door te ontwikkelen. Wij behouden ons dan ook het recht om, zonder voorafgaande kennisgeving, wijzigingen aan te brengen in deze gebruikershandleiding.